FUJIFILM Healthcare

# ELUXEO ULTRA Endoskopie System

Andreas Schuster (Clinical Marketing), Version 4.0

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


<sup>1</sup> Cancer in Germany, Robert Koch Institut, 2017

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Zahlen von IARC (International Agency for Research on Cancer)

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### Adenomas Are Detected More Often in Morning Than in Afternoon Colonoscopy

Madhusudhan R Sanaka<sup>1</sup>, Fnu Deepinder, Prashanthi N Thota, Rocio Lopez, Carol A Burke

Affiliations [+ expand](#)

PMID: 19491841 DOI: 10.1038/ajg.2009.249

[Item in Clipboard](#)

#### Abstract

**Objectives:** The effectiveness of colonoscopy in preventing colon cancer depends on adenoma detection and removal. Adequacy of bowel preparation, careful mucosal visualization, and adequate withdrawal time are known to affect adenoma detection rate (ADR). Physician fatigue, which usually increases as the day progresses, might impair ADR. The aim of this study is to assess the effect of timing of colonoscopy, morning vs. afternoon, on ADR.

**Methods:** Medical records of 9,063 colonoscopies performed in 2006 were reviewed for patient demographics, indications, timing, and findings of colonoscopy. Asymptomatic outpatients who had adequate bowel preparation and complete colonoscopy were included. Morning colonoscopies were defined as those that started before 12 noon and afternoon colonoscopies as those that started after 12 noon. ADR is defined as the detection of at least one adenoma per colonoscopy.

**Results:** A total of 3,619 colonoscopies were included, of which 1,748 (48.3%) were done in the morning and 1,871 (51.7%) were done in the afternoon. ADR was 29.3% in the morning group compared with 25.3% in the afternoon group ( $P=0.008$ ). There was a trend toward declining ADR for each subsequent hour of the day ( $P=0.01$ ). In multivariable analysis, colonoscopy in the morning was significantly associated with increased ADR (odds ratio (OR) 1.2 (1.06, 1.4)  $P=0.006$ ).

**Conclusions:** Time of performance of colonoscopy seems to be an independent predictor for adenoma detection. ADR was significantly higher in morning colonoscopies than in afternoon colonoscopies. The reasons and implications of this finding should be studied further.

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



LCI MODE

BLI MODE

WELCOME TO THE  
**BLI-PORTAL**  
YOUR ACCESS TO SCIENTIFIC INFORMATION

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For healthcare providers

## Linked color imaging reduces the miss rate of neoplastic lesions in the right colon

Editorial Supervisor: **Silvia Paggi MD**  
**Franco Radaelli MD**  
Valduce Hospital Como, Italy

**New clinical evidence**

Linked color imaging (LCI) could reduce the miss rate of neoplastic lesions in the right colon<sup>(1)</sup>.

White Light Image
LCI

**Study Center**

Valduce Hospital, Como, Italy

**Background and aim**

**Background**  
The right colon may harbor a greater prevalence of neoplastic lesions that are macroscopically more difficult to distinguish from normal mucosa, such as flat or non-polypoid lesions, that may more likely contain invasive carcinoma. Similarly, sessile serrated adenomas/polyps (SSA/P), with subtle appearance and unremarkable color, are usually located in the right colon<sup>(2)</sup>. The folds are quite deep in the right colon, and polyps are frequently located on the back of folds and on the inner curvatures of flexures, eluding the view during scope withdrawal<sup>(3)</sup>. In view of these challenges, great efforts are currently being made to develop techniques or technologies that could potentially improve polyp detection in the right colon.

**Aim**  
To investigate whether LCI improves the detection of neoplastic lesions in the right colon when compared with high definition white-light imaging (WLI).

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Figure 1



Figure 2

Type of adenoma	Endoscopic arm	Adenomas Detected at 1 <sup>st</sup> pass	Adenomas Detected at 2 <sup>nd</sup> pass	Adenoma miss rate with LCI	Adenoma miss rate with WLI
All	LCI-WLI	97	13	13 / 110 11.8%	
	WLI-LCI	77	34		34 / 111 30.6%
Advanced	LCI-WLI	6	4	4 / 10 40.0%	
	WLI-LCI	6	5		5 / 11 45.4%
SSA/Ps	LCI-WLI	4	2	2 / 6 33.3%	
	WLI-LCI	6	6		6 / 12 50.0%

[www.ncbi.nlm.nih.gov/pubmed/29539651](http://www.ncbi.nlm.nih.gov/pubmed/29539651)

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CAD e	CAD x	Title	Journal (Published)	Author
✓	✓	Performance of a new integrated CADe/CADx system for detection and characterization of colorectal neoplasia.	Endoscopy (Jan 2021)	Jochen Weigt, Alessandro Repici, Giulio Antonelli, Ahmed Afifi, Leon Kliegis, Loredana Correale, Cesare Hassan, Helmut Neumann
✓	✓	Performance of Computer-Aided Detection and Diagnosis of Colorectal Polyps Compares to That of Experienced Endoscopists	Digestive Diseases and Sciences (Aug 2021)	Taku Sakamoto, Hirotaka Nakashima, Keiko Nakamura, Ryuji Nagahama, Yutaka Saito
✓	✓	An analysis about the function of a new artificial intelligence, CAD EYE with the lesion recognition and diagnosis for colorectal polyps in clinical practice	International Journal of Colorectal Disease (Aug 2021)	Naohisa Yoshida, et al
✓		Evaluation of novel LCI CAD EYE system for real time detection of colon polyps	PLOS ONE (Aug 2021)	Helmut Neumann, et al
✓		Efficacy of a Computer Aided Detection (CADe) system in a FIT-based organized colorectal cancer screening program: a randomized controlled trial (AIFIT study).	Endoscopy (May 2022)	Rondonotti E, Di Paolo D, Rosa Rizzotto E et al.
	✓	Artificial intelligence assisted optical diagnosis for resect and discard strategy in clinical practice (Artificial intelligence BLI Characterization; ABC study)	Endoscopy (May 2022)	Rondonotti E, Hassan C, Tamanini G et al.
✓		Real-time use of artificial intelligence (CADEYE) in colorectal cancer surveillance of patients with Lynch syndrome—A randomized controlled pilot trial (CADLY)	ueg journal (Dec 2022)	Robert Hüneburg, et al
✓		Clinical Evaluation of Computer-Aided Colorectal Neoplasia Detection Using a Novel Endoscopic Artificial Intelligence: A Single-Center Randomized Controlled Trial	Digestion (Jan 2023)	Hirotaka Nakashima, et al

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

*C. Hassan, personal communication, UEGW 2019*

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	HD-WL colonoscopy	CADe-assisted colonoscopy	Risk ratio	p-value
ADR [95% CI]	45.3% [40.3-50.45%]	53.6% [48.6-58.5%]	1.180 [1.026-1.361]	—
APC (SD)	0.90 (1.32)	1.13 (1.54)	—	0.028

Endoscopist baseline (*)	ADR [95% CI]	
	HD-WL colonoscopy	CADe-assisted colonoscopy
ADR ≤ 40%	44.4% [35.3-53.9]	51.1% [42.3-59.9]
ADR 41-45%	42.8% [34.7-51.3]	53.4% [45.0-61.7]
ADR ≥ 46%	48.8% [40.0-57.7]	56.8% [47.6-65.6]

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# FUJIFILM CAD EYE + G-EYE

## First experience - 660 colonoscopies

CAD-EYE Evaluation	Standard HD Colonoscopy (Historic control)	CAD-EYE	CAD-EYE & G-EYE
Patients (n)	1362	330	330
Female (%)	64%	63%	66%
Adenoma detection rate	39,30%	<b>45%</b>	<b>63%</b>
Adenomas per patient (including serrated adenomas)	0,79	<b>1,09</b>	<b>1,31</b>
<b>Polyp Detection</b>			
Polyps (any type)		448	679
Adenomas (including serrated adenomas)		360	432
Serrated adenomas		72	130
Human eye before CAD-EYE		8%	3%
Simultaneous detection		<b>84%</b>	<b>95%</b>
Delayed human detection/interpretation		8%	2%
False positive detections (per colonoscopy)		<b>37</b>	<b>14</b>

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





All in one box

**CADEYE** Software EW10-EG01

For Desophageal Squamous Cell Carcinoma Suspect Area **LCI** **BLI**

For Gastric Neoplastic Lesion Suspect Area **White Light** **LCI** Landmark Photo Checker **White Light** **LCI**

**Expansion unit EX-1**

**CADEYE** Software EW10-EC02

CAD EYE Detection with **LCI** CAD EYE Characterisation with **BLI**

**SCALE EYE** Software EW10-EC02

FUJIFILM Holdings Corporation 19



The diagram illustrates the 'All in one box' concept for the ELUXEO system. It features a central 'Expansion unit EX-1' (a grey device) connected to various software modules. The 'CADEYE' modules include 'Software EW10-EG01' (for esophageal and gastric lesions) and 'Software EW10-EC02' (for detection and characterization). The 'SCALE EYE' module 'Software EW10-EC02' is also shown. The entire system is integrated into a mobile cart with a monitor displaying endoscopic images.

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
➤ CAD Mode Setup >> Common Setup >> CAD Mode Default >> ON/OFF


- Visual Assist Circle ON/OFF & Sound (Hi, Mid, Low, Off) can be selected.

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
## ● Movie Recording setup

➤ F9/System >> Display >> EX-WS >> Setup ↵

### ➤ **Movie Recording Setup**

- **Internal Memory: Delete latest movie file / Delete all movie files**
- **External Memory: Extract to external memory / Initialize external memory**

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